


29th August 2014

Lauren Brown

THE data processing tutorial



NEXT
GENERATION
DETECTION

Aim




- This tutorial will describe the steps required to export THE ultraFAIMS data files from Xcalibur and convert into the format required for viewing/processing data in Excel.



View data in Qual Browser


- Open Qual Browser from the Xcalibur Roadmap window.



Load RAW file in Qual Browser




- Click on the 'open data' button in the top left corner of the screen
- Select data file.



View data in Qual Browser




- TIC loads in chromatogram window as standard.



View data for ion of interest




- Ensure chromatogram window is pinned (green square indicates active window).
- Right click and select 'Ranges'.



View data for ion of interest




- Click the box under TIC.
- Change 'Plot type' to 'Mass Range' from the drop down menu.
- Enter the m/z of the ion of interest in the 'Range(s)' box
- Press 'OK'.



View data for ion of interest




- Xcalibur will only allow the export of one set of data at a time.
- By default it selects the first chromatogram range.
- TIC therefore has to be de-selected to allow the ion of interest data to be exported.
- Re-load the Chromatogram Ranges window, as the previous step, and de-select the TIC.
- Hit 'OK'.



Export data to Clipboard



- Right click on the chromatogram window and click 'Export'.
- Select 'Clipboard (Chromatogram)'. This 'copies' the data and allows it to be pasted into another program.




Export the chromatogram to the Windows clipboard

Paste data into Excel



- Owlstone will provide an Excel Spreadsheet that can be used to paste the exported data into to allow data to be viewed/processed with a compensation field axis.



Paste data into Excel



- Right click on the 'Paste data here' cell and paste in the exported data.

The screenshot displays a Microsoft Excel spreadsheet titled "T1 data processing template [Read-Only] - Microsoft Excel". The spreadsheet is divided into several columns: A, B, C, D, E, F, G, H, I, J, K, L, M, N, O, P, and Q. The data is organized as follows:

	A	B	C	D	E	F	G	H	I	J	K	L	M	N	O	P	Q
1	Ttotal	-26.6616		Time in seconds	Compensation field	CHROMATOGRAM											
2	Tcf	-2.2218				200614_ESI_1pMol_reserpine_chip002_180-280Td_-1+4Td_60s_newDFfeeder1_02.raw											
3						Data points: 921											
4	Step size	-2.25043				Time	Intensity										
5							0.00022 4895.52										
6	CF	-1		0.41239998		0.006873333	4405.703										
7				0.81130002		0.013521667	4406.505										
8				1.191		0.01985	4721.783										
9				1.59019998		0.026503333	4085.243										
10				1.97239998		0.032873333	0										
11				2.37139998		0.039523333	0										
12				2.76910002		0.046151667	0										
13				3.16690002		0.052781667	0										
14				3.56449998		0.059408333	0										
15				3.9621		0.066035	0										
16				4.35970002		0.072661667	0										
17				4.75999998		0.079333333	0										
18				5.1576		0.08596	0										
19				5.5520002		0.092586667	0										
20				5.95279998		0.099213333	0										
21				6.35059998		0.105843333	0										
22				6.74959998		0.112493333	0										
23				7.1472		0.11912	0										
24				7.54510002		0.125751667	0										
25				7.94269998		0.132378333	0										

The status bar at the bottom of the Excel window shows: Ready, FAIMS data processing template, Average: 131.8747196, Count: 1847, Sum: 242913.2335, 100%.

Convert time axis into compensation field



- The spreadsheet will automatically convert the time in minutes into time in seconds (column D).
- Manual input is then required to convert the time in seconds to compensation field (CF).
- This is due to the non-linear nature of ion-trap-type data acquisition.
- Automatic gain control means that the ion trap will fill for different amounts of time dependent on the number of ions present.
- However, the CF is ramped at a linear rate over a scan.
- This means the step in CF will not always be the same between different data points.
- However, CF across the compensation field scan can be back calculated using the steps in the following slides.

Convert time axis into compensation field



- First, the time taken for the total FAIMS scan must be calculated.
- In cell 'B1' change the equation to the cell where the FAIMS scan ends (in the example shown D914) minus the start of the FAIMS scan (D10 in the example shown).

The screenshot shows a Microsoft Excel spreadsheet titled "T1 data processing template [Read-Only] - Microsoft Excel". The formula bar at the top shows the formula for cell B1: $= (D914 - D10)$, which is circled in red. The spreadsheet data is as follows:

	A	B	C	D	E	F	G	H	I	J	K	L	M	N	O	P	Q
1	Ttotal	$= (D914 - D10)$		Time in seconds	Compensation field	CHROMATOGRAM											
2	Tcf	29.98728				200614_ESI_1pMol_reserpine_chip002_180-280Td_-1+4Td_60s_newDFfeeder1_02.raw											
3						Data points: 921											
4	Step size	0.166737				Time	Intensity										
5				0.0132			0.00022	4895.52									
6	CF	-1		0.41239998			0.006873333	4405.703									
7				0.81130002			0.013521667	4406.505									
8				1.191			0.01985	4721.783									
9				1.59019998			0.026503333	4085.243									
10				1.97239998			0.032873333	0									
11				2.37139998			0.039523333	0									
12				2.76910002			0.046151667	0									
13				3.16690002			0.052781667	0									
14				3.56449998			0.059408333	0									
15				3.9621			0.066035	0									
16				4.35970002			0.072661667	0									
17				4.75999998			0.079333333	0									
18				5.1576			0.08596	0									
19				5.5520002			0.092586667	0									
20				5.95279998			0.099213333	0									
21				6.35059998			0.105843333	0									
22				6.74959998			0.112493333	0									
23				7.1472			0.11912	0									
24				7.54510002			0.125751667	0									
25				7.94269998			0.132378333	0									

Convert time axis into compensation field



- Next, the time taken for each CF scan should be calculated.
- Divide the time taken for the FAIMS scan by the number of dispersion fields (DF) used.
 - 6 DFs were used in the example shown (cell B2).
- The number should be approximately the CF scan time inputted into the ultraFAIMS software when the data was acquired (in the example below this was 60s).

The screenshot shows an Excel spreadsheet with the following data:

	A	B	C	D	E	F	G	H	I	J	K	L	M	N	O	P
1	Ttotal	359.8474		Time in seconds	Compensation field	CHROMATOGRAM										
2	Tcf	59.97457				200614_ESI_1pMol_reserpine_chip002_180-280Td_-1+4Td_60s_newDFfeeder1_02.raw										
3						Data points: 921										
4	Step size	0.083369				Time	Intensity									
5				0.0132		0.00022	4895.52									
6	CF	-1		0.41239998		0.006873333	4405.703									
7				0.81130002		0.013521667	4406.505									
8				1.191		0.01985	4721.783									
9				1.59019998		0.026503333	4085.243									
10				1.97239998		0.032873333	0									
11				2.37139998		0.039523333	0									
12				2.76910002		0.046151667	0									
13				3.16690002		0.052781667	0									
14				3.56449998		0.059408333	0									
15				3.9621		0.066035	0									
16				4.35970002		0.072661667	0									
17				4.75999998		0.079333333	0									
18				5.1576		0.08596	0									
19				5.55520002		0.092586667	0									
20				5.95279998		0.099213333	0									
21				6.35059998		0.105843333	0									
22				6.74959998		0.112493333	0									
23				7.1472		0.11912	0									
24				7.54510002		0.125751667	0									
25				7.94269998		0.132378333	0									

Calculate compensation field step size



- Next, the time taken for each CF step should be calculated, based on the total time for the CF scan and the CF scan range.
- Divide the number of Td for the CF scan by the time taken to perform the CF scan, calculated in the previous step.
 - A 5 Td range (-1 Td to +4 Td) was used in the example shown (cell B4).

	A	B	C	D	E	F	G	H	I	J	K	L	M	N	O	P
1	Total	359.8474		Time in seconds	Compensation field	CHROMATOGRAM										
2	Tcf	59.97457				200614_ESI_1pMol_reserpine_chip002_180-280Td_-1+4Td_60s_newDFfeeder1_02.raw										
3						Data points: 921										
4	Step size	0.083369				Time	Intensity									
5				0.0132		0.00022	4895.52									
6	CF	-1		0.41239998		0.006873333	4405.703									
7				0.81130002		0.013521667	4406.505									
8				1.191		0.01985	4721.783									
9				1.59019998		0.026503333	4085.243									
10				1.97239998		0.032873333	0									
11				2.37139998		0.039523333	0									
12				2.76910002		0.046151667	0									
13				3.16690002		0.052781667	0									
14				3.56449998		0.059408333	0									
15				3.9621		0.066035	0									
16				4.35970002		0.072661667	0									
17				4.75999998		0.079333333	0									
18				5.1576		0.08596	0									
19				5.55520002		0.092586667	0									
20				5.95279998		0.099213333	0									
21				6.35059998		0.105843333	0									
22				6.74959998		0.112493333	0									
23				7.1472		0.11912	0									
24				7.54510002		0.125751667	0									
25				7.94269998		0.132378333	0									

Calculate compensation field for a given time point



- A formula is inputted into cell B6 to allow the calculation of the CF at a given time point.
- This formula uses the MOD function which tells you what is the remainder after dividing one number with another.
- This allows CF to be calculated at a particular time point based on how much time has elapsed since the start of the scan.
- The CF is calculated by taking the starting CF plus the step size (in CF s⁻¹) and multiplying it by the elapsed time.

E.g. $-1 + \$B\$4 * \text{MOD}(D10 - \$D\$10, \$B\$2)$

Calculate compensation field for a given time point



- It is important to enter the correct starting CF and step size, as well as the correct starting time. This must be changed for each exported data file.
- In the e.g. shown below, the CF scan starts at -1 Td, the step size is calculated in cell B4. The start time of the scan is in cell D10 and the total time taken is calculated in cell B2.

The screenshot shows an Excel spreadsheet titled "T1 data processing template [Read-Only] - Microsoft Excel". The formula bar for cell B6 is highlighted with a red circle, containing the formula $=-1+\$B\$4*\text{MOD}(D10-\$D\$10,\$B\$2)$. The spreadsheet data is as follows:

	A	B	C	D	E	F	G	H	I	J	K	L	M	N	O	P	Q
1	Ttotal	359.8474		Time in seconds	Compensation field	CHROMATOGRAM											
2	Tcf	59.97457				200614_ESI_1pMol_reserpine_chip002_180-280Td_-1+4Td_60s_newDFfeeder1_02.raw											
3						Data points: 921											
4	Step size	0.083369				Time	Intensity										
5							0.00022	4895.52									
6	CF	-1		0.41239998		0.006873333	4405.703										
7				0.81130002		0.013521667	4406.505										
8				1.191		0.01985	4721.783										
9				1.59019998		0.026503333	4085.243										
10				1.97239998		0.032873333	0										
11				2.37139998		0.039523333	0										
12				2.76910002		0.046151667	0										
13				3.16690002		0.052781667	0										
14				3.56449998		0.059408333	0										
15				3.9621		0.066035	0										
16				4.35970002		0.072661667	0										
17				4.75999998		0.079333333	0										
18				5.1576		0.08596	0										
19				5.55520002		0.092586667	0										
20				5.95279998		0.099213333	0										
21				6.35059998		0.105843333	0										
22				6.74959998		0.112493333	0										
23				7.1472		0.11912	0										
24				7.54510002		0.125751667	0										
25				7.94269998		0.132378333	0										

Convert each time point to compensation field



- The equation can be copied from cell B6 and pasted into the CF cell correlating to the start of the CF scan (E10 in this example).

The screenshot shows a Microsoft Excel spreadsheet titled "T1 data processing template [Read-Only] - Microsoft Excel". The formula bar at the top displays the formula: $=1+\$B\$4*\text{MOD}(D10-\$D\$10,\$B\$2)$. The spreadsheet has columns labeled A through P. Row 10 is highlighted, and cell E10 contains the value "-1", which is circled in red. The spreadsheet data is as follows:

	A	B	C	D	E	F	G	H	I	J	K	L	M	N	O	P
1	Ttotal	359.8474		Time in seconds	Compensation field	CHROMATOGRAM										
2	Tcf	59.97457				200614_ESI_1pMol_reserpine_chip002_180-280Td_-1+4Td_60s_newDFfeeder1_02.raw										
3						Data points: 921										
4	Step size	0.083369				Time		Intensity								
5							0.00022	4895.52								
6	CF	-1		0.41239998		0.006873333	4405.703									
7				0.81130002		0.013521667	4406.505									
8				1.191		0.01985	4721.783									
9				1.59019998		0.026503333	4085.243									
10				1.97239998	-1	0.032873333	0									
11				2.37139998		0.039523333	0									
12				2.76910002		0.046151667	0									
13				3.16690002		0.052781667	0									
14				3.56449998		0.059408333	0									
15				3.9621		0.066035	0									
16				4.35970002		0.072661667	0									
17				4.75999998		0.079333333	0									
18				5.1576		0.08596	0									
19				5.5520002		0.092586667	0									
20				5.95279998		0.099213333	0									
21				6.35059998		0.105843333	0									
22				6.74959998		0.112493333	0									
23				7.1472		0.11912	0									
24				7.54510002		0.125751667	0									
25				7.94269998		0.132378333	0									

Convert each time point to compensation field



- This can then be dragged to each subsequent cell – in the example below cells D10 to D914, covering the total CF scan range.
- This then gives the data needed to be able to plot CF against intensity (Column E vs. column G).

Microsoft Excel - T1 data processing template [Read-Only]


Formula bar: $=-1+\$B\$4*\text{MOD}(D914-\$D\$10,\$B\$2)$

	A	B	C	D	E	F	G	H	I	J	K	L	M	N	O	P
1	Total	359.8474		Time in seconds	Compensation field	CHROMATOGRAM										
2	Tcf	59.97457				200614_ESI_1pMol_reserpine_chip002_180-280Td_-1+4Td_60s_newDFfeeder1_02.raw										
3						Data points: 921										
4	Step size	0.083369				Time	Intensity									
5							0.00022	4895.52								
6	CF	-1		0.41239998			0.006873333	4405.703								
7				0.81130002			0.013521667	4406.505								
8				1.191			0.01985	4721.783								
9				1.59019998			0.026503333	4085.243								
10				1.97239998	-1		0.032873333	0								
11				2.37139998	-0.9667359		0.039523333	0								
12				2.76910002	-0.933580175		0.046151667	0								
13				3.16690002	-0.900416118		0.052781667	0								
14				3.56449998	-0.867268737		0.059408333	0								
15				3.9621	-0.834121351		0.066035	0								
16				4.35970002	-0.800973965		0.072661667	0								
17				4.75999998	-0.767601489		0.079333333	0								
18				5.1576	-0.734454103		0.08596	0								
19				5.55520002	-0.701306717		0.092586667	0								
20				5.95279998	-0.668159336		0.099213333	0								
21				6.35059998	-0.634995279		0.105843333	0								
22				6.74959998	-0.601731178		0.112493333	0								
23				7.1472	-0.568583792		0.11912	0								
24				7.54510002	-0.535411396		0.125751667	0								
25				7.94269998	-0.502264015		0.132378333	0								

Plot intensity against compensation field



- Plotting intensity against CF allows the CF scan data to be visualised.



Compare compensation field scans



- Different CF scans can then be plotted on the same graph to allow a comparison of their CF spectra.
- E.g.

