


All dimensions are for indication only.
Flexible hosing and loose RF wires provide flexibility on FAIMS Development Module location

REVISIONS				
ZONE	REV.	DESCRIPTION	DATE	DRAWN BY
	AAA	Original	02/10/2013	AJT

UNLESS OTHERWISE SPECIFIED ALL DIMENSIONS ARE IN mm.		 3rd Angle Projection		DO NOT SCALE DRAWING		REVISION AAA	
SURFACE FINISH: - Ra				FINISH: -		DEBUR AND BREAK SHARP EDGES	
Tolerances				TITLE:		Quick Start Development DF feeder	
Dim Precision	Tolerance	NAME	SIGNATURE	DATE		DWG NO.	01-0391
0.	+/-0.5	DRAWN	AJT	02/10/13			A4
0.0	+/-0.1	CHK'D					
0.00	+/-0.01	APPV'D					
0.000	+/-0.001	MFG					
		Q.A			MATERIAL:		
					WEIGHT:	SCALE:2:1	SHEET 1 OF 2


A

B

C

D

E




ITEM No	#	Description
A	1	Pogo Pin Locking Ring
B	1	Pogo Pin Washer
C	1	Pogo Pin PCB Assembly
D	1	Development Kit DF Feeder Base
E	1	DF Feeder Locking Nut (M16 x 1.5 Fine)
F	1	Flexible Conduit
G	1	Conduit Adaptor Male - PTFE Insert
H	1	Alternative Conduit Adaptor Male - Male
I	1	Development FAIMS Module
J	1	DF Cable Spacer
K	1	PTFE Insert
L	3	DF Feeder Wires
M	4	M3 x 12 Cap Head Screw

UNLESS OTHERWISE SPECIFIED
ALL DIMENSIONS ARE IN mm.

SURFACE FINISH: - Ra

Tolerances
Dim Precision Tolerance
0.0 +/-0.5
0.00 +/-0.1
0.000 +/-0.01
 +/-0.001



FINISH: -

DEBUR AND
BREAK SHARP
EDGES

DO NOT SCALE DRAWING REVISION

AAA

(C) This drawing is the property of Owlstone Ltd. it is to be regarded as strictly confidential and is not to be reproduced without written permission

	NAME	SIGNATURE	DATE		
DRAWN	AJT		02/10/13		
CHK'D					
APP'V'D					
MFG					
Q.A				MATERIAL:	
				WEIGHT:	

TITLE:
**Build Your Own
Development DF feeder**

DWG NO. **01-0391** A4

SCALE:2:1 SHEET 2 OF 2